

## Recombinant Human CRCP

Catalog No: #GP12837



Package Size: #GP12837-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CRCP
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-115 amino acids of human CRCP
Target Name	CGRP receptor component
Other Names	RCP; RCP9; RPC9; CGRPRCP; CGRP-RCP
Accession No.	Swissprot:O75575Gene Accession:BC040107
Uniprot	O75575
GeneID	27297;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a membrane protein that functions as part of a receptor complex for a small neuropeptide that increases intracellular cAMP levels. Alternate transcriptional splice variants, encoding different isoforms, have been characterized.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only