

Recombinant Human CREB3L2

Catalog No: #GP12838



Package Size: #GP12838-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CREB3L2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein(247FULL)
Target Name	cAMP responsive element binding protein 3 like 2
Other Names	BBF2H7
Accession No.	Swissprot:Q70SY1 Gene Accession:BC063666
Uniprot	Q70SY1
GeneID	64764;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the cAMP response element binding protein (CREB) transcription factor family. Members of this family can dimerize but form homodimers only. The encoded protein is a transcriptional activator. Translocations between this gene on chromosome 7 and the gene fused in sarcoma on chromosome 16 can be found in some tumors. Multiple transcript variants encoding different isoforms have been found for this gene.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only