

Recombinant Human CWF19L1

Catalog No: #GP12840



Package Size: #GP12840-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CWF19L1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 339-538 amino acids of human CWF19L1
Target Name	CWF19-like 1, cell cycle control (S. pombe)
Other Names	C19L1; hDm1; SCAR17
Accession No.	Swissprot:Q69YN2Gene Accession:BC008746
Uniprot	Q69YN2
GeneID	55280;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the CWF19 protein family. Mutations in this gene have been associated with autosomal recessive spinocerebellar ataxia-17 and mild mental retardation. Alternative splicing results in multiple transcript variants.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only