

Recombinant Human CYB5R4

Catalog No: #GP12842



Package Size: #GP12842-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human CYB5R4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 322-521 amino acids of human CYB5R4
Target Name	cytochrome b5 reductase 4
Other Names	NCB5OR; cb5/cb5R; dJ676J13.1
Accession No.	Swissprot:Q7L1T6Gene Accession:BC025380
Uniprot	Q7L1T6
GeneID	51167;
Storage	-20~-80°C, pH 7.6 PBS

Background

NADH-cytochrome b5 reductases are involved in desaturation and elongation of fatty acids, cholesterol biosynthesis, drug metabolism, and, in erythrocyte, methemoglobin reduction.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only