

Recombinant Human DACT3

Catalog No: #GP12857



Package Size: #GP12857-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DACT3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-167 amino acids of human DACT3
Target Name	dishevelled binding antagonist of beta catenin 3
Other Names	RRR1; DAPPER3
Accession No.	Swissprot:Q96B18Gene Accession:BC034052
Uniprot	Q96B18
GeneID	147906;
Storage	-20~-80°C, pH 7.6 PBS

Background

May be involved in regulation of intracellular signaling pathways during development. Specifically thought to play a role in canonical and/or non-canonical Wnt signaling pathways through interaction with DSH (Dishevelled) family proteins.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only