

Recombinant Human DDA1

Catalog No: #GP12862



Package Size: #GP12862-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DDA1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	DET1 and DDB1 associated 1
Other Names	PCIA1; C19orf58
Accession No.	Swissprot:Q9BW61Gene Accession:BC000615
Uniprot	Q9BW61
GeneID	79016;
Storage	-20~-80°C, pH 7.6 PBS

Background

May be involved in ubiquitination and subsequent proteasomal degradation of target proteins. Component of the DDD-E2 complexes which may provide a platform for interaction with CUL4A and WD repeat proteins.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only