

Recombinant Human DEF8

Catalog No: #GP12878



Package Size: #GP12878-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DEF8
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	differentially expressed in FDCP 8 homolog
Accession No.	Swissprot:Q6ZN54Gene Accession:BC015482
Uniprot	Q6ZN54
GeneID	54849;
Storage	-20~-80°C, pH 7.6 PBS

Background

The function of this protein remains unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only