

## Recombinant Human TAGAP

Catalog No: #GP12891



Package Size: #GP12891-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human TAGAP
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein
Target Name	T-cell activation RhoGTPase activating protein
Other Names	FKSG15; IDDM21; TAGAP1; ARHGAP47
Accession No.	Swissprot:Q8N103Gene Accession:BC015859
Uniprot	Q8N103
GeneID	117289;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the Rho GTPase-activator protein superfamily. The encoded protein may function as a Rho GTPase-activating protein. Alterations in this gene may be associated with several diseases, including rheumatoid arthritis, celiac disease, and multiple sclerosis. Alternate splicing results in multiple transcript variants encoding distinct isoforms.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only