

## Recombinant Human NUBP2

Catalog No: #GP12896



Package Size: #GP12896-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human NUBP2
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein(1-271)
Target Name	nucleotide binding protein 2
Other Names	CFD1; NBP 2; NUBP1
Accession No.	Swissprot:Q9Y5Y2Gene Accession:BC008005
Uniprot	Q9Y5Y2
GeneID	10101;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes an adenosine triphosphate (ATP) and metal-binding protein that is required for the assembly of cytosolic iron-sulfur proteins. The encoded protein functions in a heterotetramer with nucleotide-binding protein 1 (NUBP1). Alternative splicing results in multiple transcript variants.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only