

Recombinant Human ERP27

Catalog No: #GP12898



Package Size: #GP12898-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human ERP27
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 26-273 amino acids of human ERP27
Target Name	endoplasmic reticulum protein 27
Other Names	PDIA8; C12orf46
Accession No.	Swissprot:Q96DN0Gene Accession:BC030218
Uniprot	Q96DN0
GeneID	121506;
Storage	-20~-80°C, pH 7.6 PBS

Background

ERP27 is a noncatalytic member of the protein disulfide isomerase (PDI; see MIM 608012) family of endoplasmic reticulum (ER) proteins (Alanen et al., 2006 [PubMed 16940051]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only