

# Recombinant Human C11orf49

Catalog No: #GP12899



Package Size: #GP12899-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human C11orf49
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-274 amino acids of human C11orf49
Target Name	chromosome 11 open reading frame 49
Accession No.	Swissprot:Q9H6J7Gene Accession:BC009455
Uniprot	Q9H6J7
GeneID	79096;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The function of this protein remains unknown.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only