

Recombinant Human EXOSC3

Catalog No: #GP12900



Package Size: #GP12900-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human EXOSC3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-275 amino acids of human EXOSC3
Target Name	exosome component 3
Other Names	p10; PCH1B; RRP40; Rrp40p; CGI-102; hRrp-40; bA3J10.7
Accession No.	Swissprot:Q9NQT5Gene Accession:BC008880
Uniprot	Q9NQT5
GeneID	51010;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a non-catalytic component of the human exosome, a complex with 3'-5' exoribonuclease activity that plays a role in numerous RNA processing and degradation activities. Related pseudogenes of this gene are found on chromosome 19 and 21. Alternatively spliced transcript variants encoding different isoforms have been described.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only