

Recombinant Human RSU1

Catalog No: #GP12901



Package Size: #GP12901-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RSU1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-277AA amino acids of human RSU1
Target Name	Ras suppressor protein 1
Other Names	RSP-1
Accession No.	Swissprot:Q15404Gene Accession:BC005993
Uniprot	Q15404
GeneID	6251;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that is involved in the Ras signal transduction pathway, growth inhibition, and nerve-growth factor induced differentiation processes, as determined in mouse and human cell line studies. In mouse, the encoded protein was initially isolated based on its ability to inhibit v-Ras transformation. Multiple alternatively spliced transcript variants for this gene have been reported; one of these variants was found only in glioma tumors.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only