

Recombinant Human RBM11

Catalog No: #GP12905



Package Size: #GP12905-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RBM11
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 115-281 amino acids of human RBM11
Target Name	RNA binding motif protein 11
Accession No.	Swissprot:P57052Gene Accession:BC030196
Uniprot	P57052
GeneID	54033;
Storage	-20~-80°C, pH 7.6 PBS

Background

Tissue-specific splicing factor with potential implication in the regulation of alternative splicing during neuron and germ cell differentiation. Antagonizes SRSF1-mediated BCL-X splicing. May affect the choice of alternative 5' splice sites by binding to specific sequences in exons and antagonizing the SR protein SRSF1.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only