

Recombinant Human NOSIP

Catalog No: #GP12910



Package Size: #GP12910-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NOSIP
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein(1-301)
Target Name	nitric oxide synthase interacting protein
Other Names	CGI-25
Accession No.	Swissprot:Q9Y314Gene Accession:BC009299
Uniprot	Q9Y314
GeneID	51070;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene may modulate the activity and localization of nitric oxide synthase (endothelial and neuronal) and thus nitric oxide production. Alternative splicing results in multiple transcript variants that encode the same protein.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only