

# Recombinant Human DNASE1L1

Catalog No: #GP12914



Package Size: #GP12914-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human DNASE1L1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 19-302 amino acids of human DNASE1L1
Target Name	deoxyribonuclease 1 like 1
Other Names	XIB; G4.8; DNL1L; DNASEX; DNAS1L1
Accession No.	Swissprot:P49184Gene Accession:BC028092
Uniprot	P49184
GeneID	1774;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a deoxyribonuclease protein that shows high sequence similarity to DNase I. The encoded protein is localized to the endoplasmic reticulum and modified by N-linked glycosylation. Alternate transcriptional splice variants encoding the same protein have been observed.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only