

# Recombinant Human GPN2

Catalog No: #GP12917



Package Size: #GP12917-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human GPN2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-310 amino acids of human GPN2
Target Name	GPN-loop GTPase 2
Other Names	ATPBD1B
Accession No.	Swissprot:Q9H9Y4Gene Accession:BC008634
Uniprot	Q9H9Y4
GeneID	54707;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Small GTPase required for proper localization of RNA polymerase II and III (RNAPII and RNAPIII). May act at an RNAP assembly step prior to nuclear import.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only