

Recombinant Human UBXN1

Catalog No: #GP12918



Package Size: #GP12918-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human UBXN1
Brief Description	Recombinant Protein
Immunogen Description	Full length fusion protein(1-312)
Target Name	UBX domain protein 1
Other Names	2B28; SAKS1; UBXD10
Accession No.	Swissprot:Q04323Gene Accession:BC001372
Uniprot	Q04323
GeneID	51035;
Storage	-20~-80°C, pH 7.6 PBS

Background

Ubiquitin-binding protein that interacts with the BRCA1-BARD1 heterodimer, and regulates its activity. Specifically binds 'Lys-6'-linked polyubiquitin chains. Interaction with autoubiquitinated BRCA1, leads to inhibit the E3 ubiquitin-protein ligase activity of the BRCA1-BARD1 heterodimer. Component of a complex required to couple deglycosylation and proteasome-mediated degradation of misfolded proteins in the endoplasmic reticulum that are retrotranslocated in the cytosol.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only