

Recombinant Human DHRS1

Catalog No: #GP12919



Package Size: #GP12919-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human DHRS1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 2-313 amino acids of human DHRS1
Target Name	dehydrogenase/reductase 1
Other Names	SDR19C1
Accession No.	Swissprot:Q96LJ7Gene Accession:BC015943
Uniprot	Q96LJ7
GeneID	115817;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a member of the short-chain dehydrogenases/reductases (SDR) family. The encoded enzyme contains a conserved catalytic domain and likely functions as an oxidoreductase. Multiple alternatively spliced variants, encoding the same protein, have been identified.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only