

Recombinant Human TMEM248

Catalog No: #GP12921



Package Size: #GP12921-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TMEM248
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-230 amino acids of human TMEM248
Target Name	transmembrane protein 248
Other Names	C7orf42
Accession No.	Swissprot:Q9NWD8Gene Accession:BC008675
Uniprot	Q9NWD8
GeneID	55069;
Storage	-20~-80°C, pH 7.6 PBS

Background

The function of this protein remains unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only