

Recombinant Human FYTTD1

Catalog No: #GP12923



Package Size: #GP12923-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human FYTTD1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-295 amino acids of human FYTTD1
Target Name	forty-two-three domain containing 1
Other Names	UIF
Accession No.	Swissprot:Q96QD9Gene Accession:BC035006
Uniprot	Q96QD9
GeneID	84248;
Storage	-20~-80°C, pH 7.6 PBS

Background

Required for mRNA export from the nucleus to the cytoplasm. Acts as an adapter that uses the DDX39B/UAP56-NFX1 pathway to ensure efficient mRNA export and delivering to the nuclear pore. Associates with spliced and unspliced mRNAs simultaneously with ALYREF/THOC4.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only