

Recombinant Human SLC22A24

Catalog No: #GP12925



Package Size: #GP12925-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SLC22A24
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-177 amino acids of human SLC22A24
Target Name	solute carrier family 22 member 24
Other Names	NET46
Accession No.	Swissprot:Q8N4F4Gene Accession:BC034394
Uniprot	Q8N4F4
GeneID	283238;
Storage	-20~-80°C, pH 7.6 PBS

Background

SLC22A24 belongs to a large family of transmembrane proteins that function as uniporters, symporters, and antiporters to transport organic ions across cell membranes (Jacobsson et al., 2007 [PubMed 17714910]).

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only