

Recombinant Human RADIL

Catalog No: #GP12947



Package Size: #GP12947-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RADIL
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 876-1075 amino acids of human RADIL
Target Name	Rap associating with DIL domain
Other Names	RASIP2
Accession No.	Swissprot:Q96JH8Gene Accession:BC004919
Uniprot	Q96JH8
GeneID	55698;
Storage	-20~-80°C, pH 7.6 PBS

Background

Downstream effector of Rap required for cell adhesion and migration of neural crest precursors during development.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only