

Recombinant Human PISD

Catalog No: #GP12948



Package Size: #GP12948-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human PISD
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 210-409 amino acids of human PISD
Target Name	phosphatidylserine decarboxylase
Other Names	PSD; PSDC; PSSC; DJ858B16; dJ858B16.2
Accession No.	Swissprot:Q9UG56Gene Accession:BC001482
Uniprot	Q9UG56
GeneID	23761;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene catalyzes the conversion of phosphatidylserine to phosphatidylethanolamine in the inner mitochondrial membrane. The encoded protein is active in phospholipid metabolism and interorganelle trafficking of phosphatidylserine.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only