

# Recombinant Human ERGIC3

Catalog No: #GP12950



Package Size: #GP12950-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human ERGIC3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 184-383 amino acids of human ERGIC3
Target Name	ERGIC and golgi 3
Other Names	Erv46; CGI-54; C2orf47; PRO0989; C20orf47; NY-BR-84; SDBCAG84; dJ477O4.2
Accession No.	Swissprot:Q9Y282Gene Accession:BC014014
Uniprot	Q9Y282
GeneID	51614;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Possible role in transport between endoplasmic reticulum and Golgi.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only