

Recombinant Human TUFT1

Catalog No: #GP12957



Package Size: #GP12957-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human TUFT1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 191-390 amino acids of human TUFT1
Target Name	tuftelin 1
Accession No.	Swissprot:Q9NNX1Gene Accession:BC008301
Uniprot	Q9NNX1
GeneID	7286;
Storage	-20~-80°C, pH 7.6 PBS

Background

Tuftelin is an acidic protein that is thought to play a role in dental enamel mineralization and is implicated in caries susceptibility. It is also thought to be involved with adaptation to hypoxia, mesenchymal stem cell function, and neurotrophin nerve growth factor mediated neuronal differentiation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only