

Recombinant Human BCAS3

Catalog No: #GP12958



Package Size: #GP12958-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BCAS3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 729-928 amino acids of human BCAS3
Target Name	BCAS3, microtubule associated cell migration factor
Other Names	MAAB; GAOB1
Accession No.	Swissprot:Q9H6U6Gene Accession:BC001250
Uniprot	Q9H6U6
GeneID	54828;
Storage	-20~-80°C, pH 7.6 PBS

Background

Plays a role in angiogenesis. Participates in the regulation of cell polarity and directional endothelial cell migration by mediating both the activation and recruitment of CDC42 and the reorganization of the actin cytoskeleton at the cell leading edge. Promotes filipodia formation (By similarity). Functions synergistically with PELP1 as a transcriptional coactivator of estrogen receptor-responsive genes. Stimulates histone acetyltransferase activity. Binds to chromatin.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only