

Recombinant Human RMND5B

Catalog No: #GP12962



Package Size: #GP12962-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human RMND5B
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 194-393 amino acids of human RMND5B
Target Name	required for meiotic nuclear division 5 homolog B
Other Names	GID2; GID2B
Accession No.	Swissprot:Q96G75Gene Accession:BC009911
Uniprot	Q96G75
GeneID	64777;
Storage	-20~-80°C, pH 7.6 PBS

Background

The function of this protein remains unknown.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only