

Recombinant Human INPP1

Catalog No: #GP12965



Package Size: #GP12965-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human INPP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 200-399 amino acids of human INPP1
Target Name	inositol polyphosphate-1-phosphatase
Accession No.	Swissprot:P49441Gene Accession:BC015496
Uniprot	P49441
GeneID	3628;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes the enzyme inositol polyphosphate-1-phosphatase, one of the enzymes involved in phosphatidylinositol signaling pathways. This enzyme removes the phosphate group at position 1 of the inositol ring from the polyphosphates inositol 1,4-bisphosphate and inositol 1,3,4-trisphosphate.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only