

Recombinant Human OLFML3

Catalog No: #GP12970



Package Size: #GP12970-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human OLFML3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 207-406 amino acids of human OLFML3
Target Name	olfactomedin like 3
Other Names	OLF44; HNOEL-iso
Accession No.	Swissprot:Q9NRN5Gene Accession:BC009920
Uniprot	Q9NRN5
GeneID	56944;
Storage	-20~-80°C, pH 7.6 PBS

Background

Secreted scaffold protein that plays an essential role in dorsoventral patterning during early development. Stabilizes axial formation by restricting chordin (CHRD) activity on the dorsal side. Acts by facilitating the association between the tolloid proteases and their substrate chordin (CHRD), leading to enhance chordin (CHRD) degradation (By similarity). May have matrix-related function involved in placental and embryonic development, or play a similar role in other physiological processes.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only