

Recombinant Human NOLC1

Catalog No: #GP12975



Package Size: #GP12975-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NOLC1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 500-699 amino acids of human NOLC1
Target Name	nucleolar and coiled-body phosphoprotein 1
Other Names	P130; NOPP130; NOPP140; NS5ATP13
Accession No.	Swissprot:Q14978Gene Accession:BC006769
Uniprot	Q14978
GeneID	9221;
Storage	-20~-80°C, pH 7.6 PBS

Background

Related to nucleologenesis, may play a role in the maintenance of the fundamental structure of the fibrillar center and dense fibrillar component in the nucleolus. It has intrinsic GTPase and ATPase activities. May play an important role in transcription catalyzed by RNA polymerase I.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only