

Recombinant Human BZW2

Catalog No: #GP12976



Package Size: #GP12976-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	Recombinant Human BZW2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 220-419 amino acids of human BZW2
Target Name	basic leucine zipper and W2 domains 2
Other Names	MST017; HSPC028; MSTP017
Accession No.	Swissprot:Q9Y6E2Gene Accession:BC003056
Uniprot	Q9Y6E2
GeneID	28969;
Storage	-20~-80°C, pH 7.6 PBS

Background

May be involved in neuronal differentiation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only