

# Recombinant Human FAH

Catalog No: #GP12977



Package Size: #GP12977-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human FAH
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 220-419 amino acids of human FAH
Target Name	fumarylacetoacetate hydrolase
Accession No.	Swissprot:P16930Gene Accession:BC002527
Uniprot	P16930
GeneID	2184;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes the last enzyme in the tyrosine catabolism pathway. FAH deficiency is associated with Type 1 hereditary tyrosinemia (HT).?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only