

# Recombinant Human GCNT3

Catalog No: #GP12982



Package Size: #GP12982-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human GCNT3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 239-438 amino acids of human GCNT3
Target Name	glucosaminyl (N-acetyl) transferase 3, mucin type
Other Names	GNTM; C24GNT; C2GNT2; C2GNTM; C2/4GnT
Accession No.	Swissprot:O95395Gene Accession:BC017032
Uniprot	O95395
GeneID	9245;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a member of the N-acetylglucosaminyltransferase family. The encoded protein is a beta-6-N-acetylglucosamine-transferase that catalyzes the formation of core 2 and core 4 O-glycans on mucin-type glycoproteins.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only