

# Recombinant Human GRWD1

Catalog No: #GP12987



Package Size: #GP12987-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human GRWD1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 247-446 amino acids of human GRWD1
Target Name	glutamate rich WD repeat containing 1
Other Names	CDW4; GRWD; RRB1; WDR28
Accession No.	Swissprot:Q9BQ67Gene Accession:BC002440
Uniprot	Q9BQ67
GeneID	83743;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This gene encodes a glutamate-rich protein that contains five WD-repeat motifs. The encoded protein may play a critical role in ribosome biogenesis and may also play a role in histone methylation through interactions with CUL4-DDB1 ubiquitin E3 ligase.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only