

## Recombinant Human RRP1

Catalog No: #GP12996



Package Size: #GP12996-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human RRP1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 262-461 amino acids of human RRP1
Target Name	ribosomal RNA processing 1
Other Names	NNP-1; NOP52; RRP1A; D21S2056E
Accession No.	Swissprot:P56182Gene Accession:BC000380
Uniprot	P56182
GeneID	8568;
Storage	-20~-80°C, pH 7.6 PBS

## Background

The protein encoded by this gene is the putative homolog of the yeast ribosomal RNA processing protein RRP1. The encoded protein is involved in the late stages of nucleogenesis at the end of mitosis, and may be required for the generation of 28S rRNA.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only