

Recombinant Human MLLT6

Catalog No: #GP12997



Package Size: #GP12997-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human MLLT6
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 894-1093 amino acids of human MLLT6
Target Name	MLLT6, PHD finger domain containing
Other Names	AF17
Accession No.	Swissprot:P55198Gene Accession:BC007237
Uniprot	P55198
GeneID	4302;
Storage	-20~-80°C, pH 7.6 PBS

Background

AF17 was discovered as part of a genetic rearrangement in the MLL (mixed lineage leukemia) gene associated with acute myeloid leukemia. AF17 is also known as MLLT6.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only