

Recombinant Human ZBTB43

Catalog No: #GP13001



Package Size: #GP13001-1 100ug

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

| | |
|-----------------------|---|
| Product Name | Recombinant Human ZBTB43 |
| Brief Description | Recombinant Protein |
| Immunogen Description | Fusion protein corresponding to a region derived from 268-467 amino acids of human ZBTB43 |
| Target Name | zinc finger and BTB domain containing 43 |
| Other Names | ZNF-X; ZBTB22B; ZNF297B |
| Accession No. | Swissprot:O43298Gene Accession:BC008828 |
| Uniprot | O43298 |
| GeneID | 23099; |
| Storage | -20~-80°C, pH 7.6 PBS |

Background

May be involved in transcriptional regulation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only