

# Recombinant Human TBC1D2

Catalog No: #GP13002



Package Size: #GP13002-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human TBC1D2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 729-928 amino acids of human TBC1D2
Target Name	TBC1 domain family member 2
Other Names	PARIS1; PARIS-1; TBC1D2A
Accession No.	Swissprot:Q9BYX2Gene Accession:BC028918
Uniprot	Q9BYX2
GeneID	55357;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Acts as GTPase-activating protein for RAB7A. Signal effector acting as a linker between RAC1 and RAB7A, leading to RAB7A inactivation and subsequent inhibition of cadherin degradation and reduced cell-cell adhesion.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only