

Recombinant Human SAMM50

Catalog No: #GP13003



Package Size: #GP13003-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human SAMM50
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 270-469 amino acids of human SAMM50
Target Name	SAMM50 sorting and assembly machinery component
Other Names	OMP85; SAM50; TOB55; TRG-3; CGI-51; YNL026W
Accession No.	Swissprot:Q9Y512Gene Accession:BC011681
Uniprot	Q9Y512
GeneID	25813;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a component of the Sorting and Assembly Machinery (SAM) of the mitochondrial outer membrane. The Sam complex functions in the assembly of beta-barrel proteins into the outer mitochondrial membrane.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only