

Recombinant Human EIF2S3

Catalog No: #GP13005



Package Size: #GP13005-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human EIF2S3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 273-472 amino acids of human EIF2S3
Target Name	eukaryotic translation initiation factor 2 subunit gamma
Other Names	EIF2; EIF2G; MRXSBRK; eIF-2gA; EIF2gamma
Accession No.	Swissprot:P41091Gene Accession:BC019906
Uniprot	P41091
GeneID	1968;
Storage	-20~-80°C, pH 7.6 PBS

Background

The protein encoded by this gene is the largest subunit of a heterotrimeric GTP-binding protein involved in the recruitment of methionyl-tRNA(i) to the 40 S ribosomal subunit.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only