

Recombinant Human GTF2F1

Catalog No: #GP13021



Package Size: #GP13021-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human GTF2F1
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 318-517 amino acids of human GTF2F1
Target Name	general transcription factor IIF subunit 1
Other Names	BTF4; RAP74; TF2F1; TFIIF
Accession No.	Swissprot:P35269Gene Accession:BC013007
Uniprot	P35269
GeneID	2962;
Storage	-20~-80°C, pH 7.6 PBS

Background

TFIIF is a general transcription initiation factor that binds to RNA polymerase II and helps to recruit it to the initiation complex in collaboration with TFIIB. It promotes transcription elongation.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only