

Recombinant Human M1AP

Catalog No: #GP13023



Package Size: #GP13023-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human M1AP
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 380-525 amino acids of human M1AP
Target Name	meiosis 1 associated protein
Other Names	D6Mm5e; C2orf65; SPATA37
Accession No.	Swissprot:Q8TC57Gene Accession:BC025997
Uniprot	Q8TC57
GeneID	130951;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes a protein that is likely to function in progression of meiosis. A similar protein in mouse plays a role in gametogenesis in both sexes. Alternate splicing results in multiple transcript variants encoding different isoforms.?

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only