

## Recombinant Human FTSJ3

Catalog No: #GP13026



Package Size: #GP13026-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human FTSJ3
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 648-847 amino acids of human FTSJ3
Target Name	FtsJ homolog 3
Other Names	SPB1; EPCS3
Accession No.	Swissprot:Q8IY81Gene Accession:BC000131
Uniprot	Q8IY81
GeneID	117246;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Although the function of this gene is not known, the existence of this gene is supported by mRNA and EST data. A possible function of the encoded protein can be inferred from amino acid sequence similarity to the E.coli FtsJ protein and to a mouse protein possibly involved in embryogenesis.?

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only