

Recombinant Human NCOA4

Catalog No: #GP13027



Package Size: #GP13027-1 100ug

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Recombinant Human NCOA4
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 1-200 amino acids of human NCOA4
Target Name	nuclear receptor coactivator 4
Other Names	RFG; ELE1; PTC3; ARA70
Accession No.	Swissprot:Q13772Gene Accession:BC012736
Uniprot	Q13772
GeneID	8031;
Storage	-20~-80°C, pH 7.6 PBS

Background

This gene encodes an androgen receptor coactivator. The encoded protein interacts with the androgen receptor in a ligand-dependent manner to enhance its transcriptional activity. Chromosomal translocations between this gene and the ret tyrosine kinase gene, also located on chromosome 10, have been associated with papillary thyroid carcinoma. Alternatively spliced transcript variants have been described. Pseudogenes are present on chromosomes 4, 5, 10, and 14.

References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only