

# Recombinant Human AOA

Catalog No: #GP13028



Package Size: #GP13028-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human AOA
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 376-575 amino acids of human AOA
Target Name	acyloxyacyl hydrolase
Accession No.	Swissprot:P28039Gene Accession:BC025698
Uniprot	P28039
GeneID	313;
Storage	-20~-80°C, pH 7.6 PBS

## Background

This locus encodes both the light and heavy subunits of acyloxyacyl hydrolase. The encoded enzyme catalyzes the hydrolysis of acyloxylacyl-linked fatty acyl chains from bacterial lipopolysaccharides, effectively detoxifying these molecules. The encoded protein may play a role in modulating host inflammatory response to gram-negative bacteria. Alternatively spliced transcript variants have been described.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only