

STYK1 Antibody

Catalog No: #43150

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	STYK1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	WB
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total STYK1 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human STYK1
Target Name	STYK1
Other Names	NOK; SuRTK106
Accession No.	Swiss-Prot#: Q6J9G0Gene ID: 55359
Uniprot	Q6J9G0
GeneID	55359;
Calculated MW	48kd
Concentration	3.7mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:500-1:2000

Immunohistochemistry: 1:25-1:100

Images



Gel: 8%SDS-PAGE

Lysate: 40 µg

Lane: A431 cells

Primary antibody: 1/500 dilution

Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution

Exposure time: 2 minutes

Background

Receptor protein tyrosine kinases, like STYK1, play important roles in diverse cellular and developmental processes, such as cell proliferation,

differentiation, and survival.

Note: This product is for in vitro research use only