

## Recombinant Human SCM1

Catalog No: #GP13036



Package Size: #GP13036-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |   |
|-----------------------|---|
| Product Name          | Recombinant Human SCM1  |
| Brief Description     | Recombinant Protein   |
| Immunogen Description | Fusion protein corresponding to a region derived from 461-660 amino acids of human SCM1 |
| Target Name           | Scm polycomb group protein homolog 1  |
| Other Names           | Scml3   |
| Accession No.         | Swissprot:Q96GD3Gene Accession:BC009752   |
| Uniprot               | Q96GD3  |
| GeneID                | 22955;  |
| Storage               | -20~-80°C, pH 7.6 PBS   |

## Background

Associates with Polycomb group (PcG) multiprotein complexes; the complex class is required to maintain the transcriptionally repressive state of some genes.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only