

## Recombinant Human HOOK2

Catalog No: #GP13048



Package Size: #GP13048-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human HOOK2
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 518-717 amino acids of human HOOK2
Target Name	hook microtubule tethering protein 2
Other Names	HK2
Accession No.	Swissprot:Q96ED9Gene Accession:BC012443
Uniprot	Q96ED9
GeneID	29911;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Hook proteins are cytosolic coiled-coil proteins that contain conserved N-terminal domains, which attach to microtubules, and more divergent C-terminal domains, which mediate binding to organelles. The *Drosophila* Hook protein is a component of the endocytic compartment.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only