

# Recombinant Human CNOT10

Catalog No: #GP13051



Package Size: #GP13051-1 100ug

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Recombinant Human CNOT10
Brief Description	Recombinant Protein
Immunogen Description	Fusion protein corresponding to a region derived from 545-744 amino acids of human CNOT10
Target Name	CCR4-NOT transcription complex subunit 10
Accession No.	Swissprot:Q9H9A5Gene Accession:BC002931
Uniprot	Q9H9A5
GeneID	25904;
Storage	-20~-80°C, pH 7.6 PBS

## Background

Component of the CCR4-NOT complex which is one of the major cellular mRNA deadenylases and is linked to various cellular processes including bulk mRNA degradation, miRNA-mediated repression, translational repression during translational initiation and general transcription regulation. Additional complex functions may be a consequence of its influence on mRNA expression. Is not required for association of CNOT7 to the CCR4-NOT complex.

## References

Note: For in vitro research use only, not for diagnostic or therapeutic use. This product is not a medical device.

Note: This product is for in vitro research use only